


OVERBURDEN CELL PERFORMANCE COMPARISON




The Daedalus overburden cell was compared by an independent laboratory to units offered by two other major vendors of overburden cells. The tests were conducted using an Oxford GeoSpec2 at 2 MHz or 12 MHz employing CPMG experiments used for porosity measurements. The Daedalus unit tested was a standard 49 mm OD cell for 1” diameter cores rated to 10,000 psi.

**At 2 MHz:** It is clear the Daedalus cell trace (red) has a much lower background signal than either vendor 1 (orange) or vendor 2 (blue). Only the Daedalus cell would permit measurements of shale at short TE values; the other cell signals would swamp the sample signal leading to large errors in measurements.



**At 12 MHz:** Again, it is readily apparent that the background signal from the Daedalus cell (red) is much lower than the vendor 1 cell (orange).



The cell from vendor 2 could not be tested because it was not designed to work at that field with the instruments available for this testing.




200 Racoonin Drive, Suite 106  
Aston, PA 19014  
United States


Phone: 610-358-4728  
Fax: 610-361-8509  
E-mail: sales@daedalusinnovations.com

Accurate measurements of porosity data from tight rocks (shale) requires short time-to-echo spacings. Further, the cumulative pore volume of shale samples tends to be small. Only the Daedalus cell has sufficiently low background signal to permit high quality measurements of such samples. The cells from vendor 1 & 2 require operation at longer TE values thus throwing away data available at short  $T_2$  values (small pores) thus decreasing the accuracy of the results obtained. A comparison between a shale sample in a Daedalus cell and that a cell from vendor 1 is shown on the opposite side.

**Shale measurements:** Porosity data was collected on shale samples in both the Daedalus cell and the cell from vendor 1 at 12 MHz. For comparison purposes the sample was also measured in a glass vial to assess the accuracy of the data collected. The method of data collection was the same as that used for measurement of the background signal from the empty cells.



The low background signal from the Daedalus cell allows data to be collected with  $TE = 114 \mu s$ . The green traces are the data of the sample in a glass vial with the red traces the data from the Daedalus cell. Though the Daedalus cell data does overestimate the pore volume slightly the overall agreement is remarkably good. No effort was made to subtract the empty cell spectrum from the sample spectrum in this analysis.



The data from the vendor 1 cell could not be collected with  $TE = 114 \mu s$  due to the background signal overwhelming the sample data signal. Even so, at  $TE = 200 \mu s$  the total pore volume is considerably overestimated (orange). Further, due to the high background signal at low  $T_2$  values the overall trajectory of the trace is influenced such that the overall profile only marginally resembles the glass vial data (green).



200 Racoosin Drive, Suite 106  
Aston, PA 19014  
United States

Phone: 610-358-4728  
Fax: 610-361-8509  
E-mail: sales@daedalusinnovations.com

When the background signal from core holder is too large it may inhibit measurements of shale samples, or risk throwing away valuable data at low  $T_2$  values. With the Daedalus cell such measurements can be performed with a high accuracy at short TE values. The signal is sufficiently low that it does not overwhelm the sample signal such that it makes sense to subtract the background signal from the sample to obtain even more accurate results. Such methods are currently being explored.